www.mobilehealth2u.com

LEVER CHEK[®]
Master your life.

TaiDoc Technology Corporation
B1-7F, No. 127, Wugong 2nd Rd.,
Wugu Dist., 24888 New Taipei City, Taiwan
www.taidoc.com

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LEVER CHEK[®]
TD-3250

2_{in}**1**

Blood Glucose plus
Blood Pressure Monitoring
System

OWNER'S MANUAL

**ACCURATE
SMART
PRECISE**

Blood Glucose plus Blood Pressure Monitoring System

This owner's manual contains important information that you must know about your system. Please read it carefully and keep it for future reference.

For other questions regarding this system, please contact your local customer service. At all other times, you should contact your health care professional for assistance.

IMPORTANT SAFETY INSTRUCTIONS

READ THIS BEFORE USING

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
2. Use the device only for the intended use described in this manual.
3. Do not use accessories which are not supplied by the manufacturer.
4. Do not let the equipment or its flexible cord come into contact with surfaces which are too hot to touch.
5. Do not use the equipment where aerosol sprays are being used, or where oxygen is being administered.
6. Do not use device if it is not working properly, or if it has suffered any damage.
7. Before using product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

KEEP THESE INSTRUCTIONS

Table of Contents

BEFORE YOU START	05
Warnings and Precautions	05
Health Information	06
ABOUT THIS SYSTEM	08
Intended Use	08
Principle of Measurement	08
Major Features	08
Content of the System	10
Appearance and Key Functions of the Monitor	12
LCD Display	14
Appearance of the Test Strip	15
PREPARATIONS BEFORE USE	16
Battery Replacement	16
Setting the Monitor	18
BLOOD GLUCOSE MEASUREMENT	20
Important Information	20
Checking the Code Number/ Coding the Monitor (For TD-3250B only)	21
Checking with TaiDoc Control Solutions	23
Important Control Solution Information	24
Doing a Control Solution Test	25
Testing with Blood Sample	28
Preparing the Puncture Site	29
Special Messages	33
Comparing Monitor and Laboratory Results	34
BLOOD PRESSURE MEASUREMENT	36
Suggestions before Measuring	36
Correct Application of the Pressure Cuff	36
Testing Your Blood Pressure	37

MONITOR MEMORY	39
Viewing Memory on the Monitor	39
Viewing Results on a Personal Computer	41
Deleting Memory	43
TAKE CARE OF YOUR MONITOR & STRIP	44
Cleaning	44
Storage	44
PROBLEM-SOLVING GUIDE	47
Error Message	48
Problem in Operation	49
SPECIFICATIONS	53
SYMBOLS INFORMATION	54
FCC STATEMENT	55

BEFORE YOU START


⊙ Warnings and Precautions

- ▶ The CLEVER CHEK TD-3250 System is designed for use on individuals **age 16 and above**. It shall **NOT** be used:
 - To diagnose newborns with diabetes.
 - On infants or persons who cannot communicate.

- ▶ This device does **NOT** serve as a cure of any symptoms or disease. The data measured are only for reference. Always consult your physician to have the results interpreted.

- ▶ This device is **NOT** able to measure in the presence of common arrhythmia, such as arterial or ventricular premature beats or arterial fibrillation. It may produce reading error.

- ▶ Do **NOT** use the device for purposes other than measuring blood glucose and blood pressure for human beings.
- ▶ Do **NOT** wrap the pressure cuff on anything besides your arm.

◎ Health Information

Blood Glucose

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that **keeping blood glucose levels close to normal** can reduce the risk of diabetes complications by up to 60%*1. The results you get with the **CLEVER CHEK TD-3250** system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Plasma glucose range (mg/dL) for people with diabetes	Your target range (mg/dL)
Fasting and before meal	70-130 mg/dL	
2 hours after meals	Less than 180 mg/dL	


Source: American Diabetes Association (2008). Standards of Medical Care in Diabetes. Diabetes Care, 31 (Supplement 1): S12–S54.

*1:ADA Clinical Practice Recommendations 2003

Please work with your doctor to determine a target range that works best for you.

Blood Pressure

Clinical studies show that the adult diabetes are often accompanied by elevated blood pressure. People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment*2. Knowing your routine blood pressure trend tells whether your body is in good condition or not. Human blood pressure naturally increases after reaching middle age. This symptom is a result of continuous aging of the blood vessel. Further causes include obesity, lack of exercise, and cholesterol (LDL) adhering to the blood vessels. Rising of blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction. The WHO (world health organization) published the guideline of blood pressure range:

Source: 1999 WHO/ISH guidelines for the management of hypertension
 *2: American Diabetes Association: The Diabetes-Heart Disease Link Surveying Attitudes, Knowledge and Risk (2002)

ABOUT THIS SYSTEM

⦿ Intended Use

The CLEVER CHEK TD-3250 Blood Glucose plus Blood Pressure Monitoring System is indicated for the quantitative measurement of glucose in fresh capillary whole blood for self testing by persons with diabetes in the home or by healthcare professionals in healthcare facilities. The system is also intended to be used to measure non-invasively the systolic and diastolic blood pressure and pulse rate of an adult individual, over age 16, at home.

⦿ Principle of Measurement


Blood glucose is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The monitor measures the current and displays the corresponding blood glucose level. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Blood pressure is measured non-invasively at the wrist based on the Oscillometric method.

Both functions work separately (one measurement either blood glucose or blood pressure at a time) to avoid any interference problems.

⦿ Major Features


Two subtype of TD-3250 were developed and the difference of function between them is simply in **data transmission**.

Subtype of TD-3250	Data transmission function
<p data-bbox="75 227 223 258">TD-3250B</p>  <p>The diagram shows a handheld device (TD-3250B) on the left and a desktop computer (PC) on the right. A cable is connected between the device and the PC, indicating a wired connection.</p>	<p data-bbox="576 227 917 335">You can view results on the PC through cable connection.</p>
<p data-bbox="75 579 223 609">TD-3250C</p>  <p>The diagram shows a handheld device (TD-3250C) on the left and a desktop computer (PC) on the right. The device has a signal icon above it, representing a wireless connection to the PC.</p>	<p data-bbox="576 579 938 725">You can view results on the PC through wireless connection using Bluetooth technology.</p>

Content of the System

These products have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Use only CLEVER CHEK TD-3250 test strips and TaiDoc control solution with your CLEVER CHEK TD-3250 Monitor.


① A blood glucose plus pressure monitor	⑧ First time user guide
② Pressure cuff	⑨ Owner's manual
③ 25 test strips	⑩ Carrying bag
④ Code strip (For TD-3250B only)	⑪ Warranty card
⑤ 2 levels of control solution	⑫ Bluetooth Adapter for model TD-3250C
⑥ 25 lancets	⑬ RS232 Cable for model TD-3250B (optional)
⑦ Lancing device	⑭ Daily logbook



CAUTION

1. Check your system to be sure that it is unopened prior to use and that it contains all parts listed above. If either of these conditions occurred, please return your system to the place of purchase.
2. Please note that RS232 Cable for model TD-3250B is optional. It is not included in the standard kit. Please ask your local agent for help.

Appearance and Key Functions of the Monitor

LCD Display

TEST RESULT AREA

Test results are displayed here. (A decimal point appears when the unit of measurement is millimoles per litre.)

DATE & TIME

BLOOD DROP SYMBOL

This symbol tells you when to apply the sample.

KETONE WARNING


Appears when glucose result is above 240mg/dL (13.3mmol/L).

TEST STRIP MESSAGE

When power of monitor is OFF and insert strip or power turn off automatically with test strip inserted in test slot, this message flashes.

MEMORY

Indicates a test result stored in memory mode.

CODE

Appears with the code number of the test strips currently in use.

SYSTOLIC PRESSURE SYMBOL

UNITS OF MEASUREMENT

mmol/L or mg/dL (mmHg or Kpa) will appear with the test result.

DIASTOLIC PRESSURE SYMBOL

BATTERY SYMBOL

Appears when the battery is low or must be replaced.

TEMPERATURE MESSAGE

Appears with the ambient temperature after inserting a test strip.

DAY AVERAGE
Indicates average of blood glucose reading measured.

PULSE RATE

Displays number of heart beats per minute (bpm).

PULSE SYMBOL

Flashes when measuring blood pressure.

◎ Appearance of the Test Strip

Your system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole of the test strip and is automatically drawn into the reaction chamber where the reaction takes place.

The test strip consists of the following parts:

Contact Bars

Insert this end of the test strip into the monitor.

Push it in firmly until it will go no further.

Test Strip Handle

Hold this part to insert the test strip into the slot.

Confirmation Window

This is where you confirm if enough blood

Has been applied to the absorbent hole of the strip.

Absorbent Hole

Apply a drop of blood here,

The blood will be sucked automatically.


See the section, Testing Your Blood Glucose, for complete instructions.

PREPARATIONS BEFORE USE


◎ Battery Replacement

Your monitor comes with batteries already installed. When replacing, use ONLY 1.5V AA size alkaline batteries for best performance and longest life. The monitor will remind you when the power is getting low by displaying two different messages:


1. Battery Symbol  appears while you perform tests. This means it is about time to change the batteries, but this testing result is still accurate.


2. Battery Symbol  flashes by itself and the monitor cannot perform any measurement. This means you must change the batteries immediately.

Make sure the monitor is off when replacing the batteries.

STEP 1- Press the edge of the battery cover and pull up. Lift the battery cover and remove the batteries inside.

STEP 2- Insert four 1.5 V AA size alkaline batteries. Make sure the + (positive) and - (negative) marks match as indicated in the battery compartment. If inserted correctly, you will hear a “beep”.

CAUTION

1. Do not use different type, date or brand name batteries together. Use only new batteries of the required size and type.
2. Replacing the batteries does not affect previous test results stored in the monitor. But you may need to update the settings.
3. Batteries might leak chemicals if not used for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e. 3 months or more).
4. As with all small batteries, the batteries should be kept away from small children who still put things in their mouths. If they are swallowed, promptly see a doctor for help.
5. Discard batteries according to your local regulations.

◎ Setting the Monitor

Time, date and units can be set in the monitor. Two keys are involved: **S** button and **M** button.

Step1 Press the S button to start setting.

Step2 Year flashes first.

Step3 Press the M button to make changes. Press the S button to next.

Step4 Date flashes. Follow Step3 for changes.

Step5 Time flashes. Follow **Step 3** for changes.

Time 10:08 PM

Step6 Unit of glucose concentration*³ flashes. Follow **Step 3** for changes.

mg/dL
mmol/L

Step7 Unit of blood pressure flashes. Follow **Step 3** for changes.

mmHg
KPa

Step8 Unit of Temperature flashes. Follow **Step 3** for changes.

°C
°F

NOTE

- The time, date, unit of measurement and unit of temperature can **ONLY be changed** in the setting mode. Therefore, when you perform a glucose testing or blood pressure testing, those parameters are not possible to be changed.
- *3: The milligram per deciliter (mg/dL) is the standard unit in the United States. The millimole per liter (mmol/L) is the standard unit in Canada. Use of the wrong unit of measure may cause you to misinterpret your blood glucose level, and may lead to incorrect treatment.
- While the monitor is in the setting mode, if no button is pressed for 3 minutes, the monitor will turn off automatically.
- You may also press the  button any time to turn off the monitor and exit the Setting Mode.

BLOOD GLUCOSE MEASUREMENT

⊙ Important Information

- Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If you get your blood glucose results lower or higher than usual, and do not have symptoms, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only capillary whole blood sample to test your blood glucose. Using other substances will cause wrong results.
- If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose plus blood pressure monitor.
- Please refer to your test strip package insert for additional important information.

⦿ **Checking the Code Number/ Coding the Monitor (For TD-3250B only)**


When using the system for the first time, or starting to use a new lot of test strips, you will check the code number. It is important to make sure that the LCD displayed code is the same as the code on the test strip vial before you proceed. **Only test blood glucose when the codes are identical, so the test results will be accurate.**

Calibration of TD-3250B monitor.

CLEVER CHEK TD-3250B monitor should be calibrated by inserting the code strip. You need to calibrate your monitor according to the following steps.

STEP 1- Fully insert the code strip into the test slot. A single beep will sound and the code number will appear on the display. Verify that the code number is the same as the code on the test strip vial.

Make sure the code number on LCD display, on the code strip, and on strip vial are the same.

STEP 2- Remove the code strip, the LCD display will show "OK". This tells you that the monitor has finished coding and is ready for blood glucose testing.

NOTE

The code number on the display is only for your reference, it may not be actual code for this meter.

WARNING

It is important to make sure that the LCD displayed code is the same as the code on the test strip vial before testing. Failure to do so will get inaccurate results.

⊙ **Checking with TaiDoc Control Solutions**

TaiDoc control solutions contain a known amount of glucose that reacts with test strips. By comparing your control solution test results with the expected range printed on the test strip vial label, it is able to check that the monitor and the test strips are working together as a system and that you are performing the test correctly. It is very important that you do this simple check routinely to make sure you get accurate results.

How often should the control solution test be performed?

- When you use this system to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test your blood.
- For routinely check the monitor and test strips, perform a single test for each level of control solution at least once a week.

When should the control solution test be performed?

- When you first get your monitor.
- When you begin using a new vial of test strips.
- Whenever you suspect that the monitor or test strips are not working properly.
- When your blood glucose test results are not consistent with how you feel, or when you think your results are not accurate.
- When your test strips are exposed to extreme environmental conditions (See **Storage** section of this manual).
- When you want to practice running the test.
- If you drop the monitor.

◎ Important Control Solution Information

- Use only TaiDoc control solutions.
- Check the expiry date on the control solution vial. Do not use if expired.
- Control solution, monitor, and test strips should come to room temperature (20~25°C/68~77°F) before testing.
- Shake the vial, discard the first drop of control solution, and wipe off the dispenser tip to ensure a good sample and an accurate result.
- Use only for 90 days after first opening. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.
- Store the control solution tightly closed at temperatures 2 to 30°C (36 to 86°F). Do not freeze.

NOTE

The control solution range printed on the test strip vial is for TaiDoc control solution only. It is used to test monitor and test strip performance. It is not recommended range for your blood glucose level.


⊙ Doing a Control Solution Test

TAKE A TEST STRIP OUT WITH CLEAN AND DRY HANDS FIRST.

STEP 1- Insert test strip

Insert a test strip with contact bars end first and facing up into the test slot. The monitor turns on automatically and displays the followings in sequence:

“CHK” and “” → ambient temperature → code number and flashing “”

Be sure the code number on the display is the same as the code number on the test strip vial. If the code numbers do not match, please see “Checking the Code Number/ Coding the Monitor” section of this manual.

STEP 2- Press the M button

While the “” symbol appears on the display, press the M button and then “CTL” will appear on the display. With the “CTL” sign on the display, the monitor will not store your test result in memory. If you decide not to perform a control solution test, press the M button again, and the “CTL” sign will disappear.


WARNING

1. Contact bars must be inserted all the way into the monitor or you may get an inaccurate test result.
2. Every time you perform a control solution test, you must enter into the “**CTL**” test mode so that the test result will not be stored in the monitor memory. Failure to do so will confuse the blood glucose test result with the control solution test result in memory.


STEP 3- Obtain control solution

Shake the control solution vial well. Remove the cap from the control solution bottle. Place cap on flat surface. Squeeze the vial, discard the first drop, and wipe off the dispenser tip to prevent contamination. Squeeze the vial again to get another drop and **apply the drop to the top of cap.**


**STEP 4- Apply Control Solution**

Hold the monitor to meet the absorbent hole of the test strip and the drop will be automatically drawn into the test strip. Make sure the confirmation window fills completely. The monitor begins to count down.

To avoid contaminating the control solution with the content of the test strip, you have to place a drop of control solution on a clean surface. Do not directly apply control solution into a strip.

**STEP 5- Read and Compare the Result**

After counting to 0, the test result of control solution is shown on the screen. Compare this result with the range printed on the test strip vial. It should fall within this range.


Out-of-range results

If test results fall outside the range printed on the test strip vial, check the section of “Problem in Operation” in troubleshooting guide and repeat the test. If you continue to get out-of-range results, it means that the system may not be working properly. Do NOT test your blood. Please contact your local customer service for help.

⊙ Testing with Blood Sample

Overview of the Lancing Device

If your lancing device differs from the one shown here, please refer to the manufacturer's manual to ensure proper usage.





WARNING

To reduce the chance of infection:

- Never share a lancet or the lancing device with anyone.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

STEP 1 Setting up the Lancing Device

1. Pull off the cap of the lancing device.
2. Insert a lancet into the lancet holder and push down firmly until it is fully secured.
3. Twist the protective disk off the lancet.
4. Replace the cap by aligning the arrow on the cap with the release button.
5. Select the depth of penetration by turning the adjustable tip in either direction so that the arrow on the cap points to the desired depth.
6. Pull the cocking control back until it clicks. You will see a color change inside the release button when it is ready.

**⊙ Preparing the Puncture Site**

Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the glucose value obtained.

Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.


Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site either at fingertips or another body parts (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and let it air dry.
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting up the lancing device.


STEP 2- Insert the test strip and check the code number

Insert a test strip with contact bars end first and facing up into the test slot. The monitor turns on automatically and displays the followings in sequence:

“CHK” and “” → ambient temperature
→ code number and flashing “”

Make sure the code number showed on the screen is the same as the code number printed in the test strip vial. If the code numbers are different, please refer to section of “Checking the Code Number/Coding the Monitor” for the procedure of coding.


STEP 3- Get a drop of blood

Select the puncture site in finger. Clean the puncture site with 70% alcohol cotton and **let it air-dry**.

Select a puncture site on fingertip

Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete.

After penetration, discard the first drop of blood with a clean tissue paper or cotton. Then gently squeeze the punctured area to obtain blood. But be attention **NOT to smear the blood sample**.


The volume of blood sample must be at least 0.7 microliter (• actual size).

**WARNING**

- Choose a different spot each time you test. Repeated puncture in the same spot may cause soreness and calluses.
- Since the first drop of blood usually contains tissue fluid and serum, which may affect the test result, it is recommended to be discarded.

STEP 4- Apply blood into the test strip

When “” is flashing on the screen, apply your blood to the absorbent hole of the test strip until the confirmation window is **fully covered** with blood. The monitor then begins to count down automatically.


STEP 5- Obtain an accurate result in 7 seconds

The result of your blood glucose test is shown after the monitor counts to 0. This reading is automatically stored in the monitor.

STEP 6- Remove the lancet

Remove the lancing device cap and the lancet. Place the disk on a hard surface and push the exposed needle tip into the protective disk. **Always use caution when removing the lancet.**

WARNING

The used lancet and the used test strip may be potentially bio-hazard. Please discard it carefully according to your local regulations.

Special Messages

MESSAGE	WHAT IT MEANS	ACTION
	<p>Appears when your result is below measurement limit, which is less than 20 mg/dL (1.1 mmol/L).</p>	<p>This indicates hypoglycemia (low blood glucose.)</p> <p>You should immediately treat hypoglycemia as recommended by your healthcare professional.</p>
	<p>Appears when your result is above measurement limit, which is higher than 600 mg/dL (33.3 mmol/L).</p>	<p>This indicates severe hyperglycemia (high blood glucose).</p> <p>You should seek immediate medical assistance.</p>
<p>240mg/dL</p> <p>KETONE</p>	<p>KETONE appears when your result is equal or higher than 240 mg/dL (13.3 mmol/L).</p>	<p>This indicates there is a possibility of ketone accumulation if you are Type 1 diabetes.</p> <p>Please seek immediate medical assistance.</p>

◎ Comparing Monitor and Laboratory Results

The monitor provides you with whole blood equivalent results. The result you obtain from your monitor may differ somewhat from your laboratory result due to normal variation. Monitor results can be affected by factors and conditions that do not affect laboratory results in the same way (See test strip package insert for typical accuracy and precision data and for important information on limitations). To make an accurate comparison between monitor and laboratory results, follow the guidelines below.

Before you go to the lab:

- Perform a control solution test to make sure that the monitor is working properly.
- It is best to fast for at least eight hours before doing comparison tests.
- Take your monitor with you to the lab.

While at the lab:

Make sure that the samples for both tests (the monitor test and the lab test are taken and tested within 15 minutes of each other).

- Wash your hands before obtaining a blood sample.
- Never use your monitor with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication, or experienced stress*⁴. In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test*⁵. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (severe dehydration) may also cause a monitor result to be different from a laboratory result.

References

*4: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*5: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.


BLOOD PRESSURE MEASUREMENT

⊙ Suggestions before Measuring


1. Avoid caffeine, tea, alcohol, and cigarette at least 30 minutes before measurement.
2. Wait 30 minutes after exercising or bathing before measurement.
3. Sit or lie down for at least 10 minutes before measuring.
4. Do not measure when feeling anxious or tense.
5. Take a 5-10 minute break between measurements. This break can be longer if necessary, depending on your physical conditions.
6. Keep the records for your physician as reference.
7. Blood pressure varies from two hands naturally. Always measure your blood pressure on the same arm.

⊙ Correct Application of the Pressure Cuff


STEP 1. Connect the air plug of the tubing to the air jack at the side of the monitor.

STEP 2. The smooth cloth is on the inside of the cuff loop and the metal D-ring will not touch your skin.

STEP 3. Hold your left arm in front of you with your palm facing up. Slide the cuff onto your left arm, above your elbow. The red line on the edge of the cuff should be approximately **0.8 to 1.2 inch (2 cm to 3 cm)** above your elbow. Align the tubing over the main arteries at the inside of your arm.


After the cuff pressure is reached, the cuff begins to deflate. You will see number decreasing and ♥ flashing in the display.

STEP 5 Read a result.

After the measurement, the monitor displays the systolic pressure, diastolic pressure, and heart rate. Press \odot button to switch off. Or it will switch off automatically after 3 minutes.

NOTE

- If a higher pressure value is needed, the monitor will stop deflation and inflate again.
- If you press the \odot button during the measurement, the monitor will turn off immediately.

MONITOR MEMORY


Your monitor stores the most 352 recent results, and also provides the average of the blood glucose test results currently stored. Please follow the steps listed below to review those results.

⊙ View Results on the Monitor

Step1 Enter the memory mode

When the monitor is off, press the M button. The 7-day average will appear, indicating that you are in the memory mode. If you continue to press the M button, the 14-, 21-, 28-, 60-, and 90-day average will appear. You can then review the average tests in the memory.

When using the monitor for the first time, “---” appears, showing that there are no test result in memory.


The 7-day average is calculated from the blood glucose results obtained during the last 7 days. The example displayed means that the blood glucose average performed over the last 7 days was 128 mg/dL.

Step2 Recall test results

After 90-day average, the other results will be displayed in the order of time and date taken. You will see the number of the result stored, followed by the result. When the memory is full, the oldest results are deleted as new ones are added.

Step3 Exit the memory mode

Press the \odot button to turn off the monitor and exit the memory mode.

NOTE

1. The control solution results are **NOT** stored in the memory (please follow the page 26 **WARNING** for information). The list of past results and the result average are for blood glucose results only.
2. If you'd like to exit from the memory, press the \odot button to turn off the monitor or wait without pressing for 3 minutes for auto-shut down.

◎ View Results on a Personal Computer

Results in memory can be transmitted to the personal computer by either cable or wireless connection for model TD-3250B or TD-3250C, respectively. Accessories needed to activate this function are:


- Health Care System Software: a kind of software downloaded from TaiDoc website.
- Interface Cable: an optional accessory for model TD-3250B.
- Bluetooth Adapter: an accessory for model TD-3250C.

If the above accessory is optional, please contact your local customer service.

Transmitting data via cable (TD-3250B)

STEP 1- Install Health Care System Software on your computer by following the instructions provided on TaiDoc's website: <http://www.taidoc.com>.

STEP 2- Connect to Personal Computer with interface cable
When the monitor is off, connect the serial port of PC with the Data Port of monitor by interface cable. Then "PCL" and time will appear in the monitor, indicating it is ready to transmit data.

STEP 3- Follow the instructions provided in the software to transmit the data. Results transmitted will include date and time. Remove the cable and the monitor will automatically turn off.

Transmitting data via Bluetooth adapter (TD-3250C)

STEP 1- Install Bluetooth Adapter Software.

STEP 2- Reset TD-3250C monitor.

STEP 3- Register TD-3250C.

STEP 4- Install Health Care System Software on your computer by following the instructions provided on TaiDoc's website:
<http://www.taidoc.com>.

STEP 5- Connect TD-3250C with PC and view your data. Please refer to the additional package insert "How to Start Wireless Communication of TD-3250C Blood Glucose plus Blood Pressure (BGBP) Monitor" for detail information.

NOTE

While the monitor is connected to the PC, it is unable to perform a blood glucose or blood pressure test.

⦿ Deleting Memory

To clear the memory, with the monitor turned off, press and hold the M button for 3 seconds. “CLr/ALL” then appears and indicates that **ALL** the memories stored in the monitor will be cleared.

TAKE CARE OF YOUR MONITOR & STRIP


To avoid the monitor and test strips getting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

⦿ Cleaning


- To clean the monitor exterior, wipe with a cloth moistened with tap water or a mild cleaning agent, then dry the monitor with a soft and dry cloth. Do not flush with water.
- Do not use organic solvents to clean the monitor and cuff.
- The cuff can be cleaned by wiping with a moistened cloth and soap. Do not immerse in water.

⦿ Storage


1. Monitor Storage

- Storage condition: $-20^{\circ}\text{C}\sim 60^{\circ}\text{C}$ ($-4^{\circ}\text{F}\sim 140^{\circ}\text{F}$), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.

- Avoid dropping and strong impact.
- Avoid direct sunlight and humidity.
- Do not disassemble, modify or try to repair the monitor or wrist cuff by yourself.


- Do not over-wring the cuff and turn the cuff inside out.
- If you are not going to use the monitor for an extended period, please remove the batteries.

2.Strip Storage


- Storage condition: 4°C~40°C (39.2°F~104°F), below 85% relative humidity. **Do NOT freeze.**
- Store your test strips in their original vial only. Do not transfer to other container.

- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately replace the vial cap and close it tightly.

- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the discard date (the date opened plus 90 days) on the vial label when you first open it. Discard remaining test strips 90 days after first opening date.

- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

3. Control Solution Storage

- Storage condition: Store the control solution tightly closed at temperatures 2 to 30°C(36 to 86°F). **Do NOT freeze.**
- Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.

PROBLEM-SOLVING GUIDE

Following is a summary of some display messages. These messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message. In the event of a problem, refer to the information under action.

Never try to disassemble the monitor in any circumstances. If you encounter any error messages not listed below or if you have followed the actions recommended below but the problem keeps unsolved, please call the customer service for support.

◎ Error Message

MESSAGE	CAUSE	WHAT TO DO
Err.00	Weak pulse.	Refit cuff tightly, relax, and repeat measurement as shown in "Correct Application of the Pressure Cuff".
Err.02	Monitor can not figure out systolic pressure.	
Err.04	Monitor can not figure out diastolic pressure or it is out of range.	
Err.22	Invalid coding.	Try a new code strip and insert again.
Err.24	Used strip insertion.	Use a new strip.
Err.25	Ambient temperature is below operating temperature.	Please operate this monitor between 10°C ~40°C(50°F~104°F)
Err.26	Ambient temperature is below operating temperature.	
Err.05 Err.06 Err.09 Err.10 Err.11 Err.15 Err.20 Err.28	<p>Those messages indicate internal problems with the monitor. You cannot fix by yourself.</p> <p>Please contact Customer Service Line for help.</p>	

⊙ Problem in Operation

1. Blood Glucose Measurement

What happened?

The monitor does not display a message after inserting a test strip.

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.
Test strip inserted upside down or not completely inserted.	Insert the test strip correctly with the contact bars end first and facing up.
Defective monitor.	Contact local customer service.

What happened?

The test does not start after applying the sample.

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample.	Replace the batteries. Repeat the test with a new test strip and a larger sample.
Defective test strip.	Repeat the test with a new test strip.
Sample was applied at the time that not flashing ◀ appears on the display.	Repeat the test with a new test strip. Apply sample only when ◀ appears on the display.
Defective monitor.	Contact local customer service.

What happened?

If you turned off the monitor but it makes a beeping sound.

POSSIBLE CAUSE	WHAT TO DO
The strip is still inside the slot.	Remove the strip. If you need to test blood glucose, insert an unused strip.

What happened?

If the control solution test result is out of range.

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read the instruction thoroughly and repeat the test again.
Improper code number.	Check if the code number on the display matches the code number on the test strip vial.
Do not shake the control solution vial very well.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date or the discarded date of the control solution.
Control solution that is too warm or too cold.	Control solution, monitor, and test strips should come to room temperature (20~25°C/68~77°F) before testing.
Test strip deterioration.	Repeat the test with a new test strip.
Monitor malfunction.	Contact local customer service.

2. Blood Pressure Measurement

What happened?

No display after pushing the  button.

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

What happened?

Heart rate is higher/lower than user's average.

POSSIBLE CAUSE	WHAT TO DO
Moving during measurement.	Repeat measurement.
Measuring right after exercise.	Rest at least 30 minutes before measurement.

What happened?

A result is higher/lower than user's average measurement.

POSSIBLE CAUSE	WHAT TO DO
May be not in correct position during measuring.	Adjust to a correct position to measure.
Blood pressure naturally varies from time to time.	Keep in mind for next measurement.

What happened?

Cuff inflates again during measuring.

POSSIBLE CAUSE	WHAT TO DO
Cuff is not fastened.	Fasten the cuff again.
Normal action. If user's blood pressure is higher than the previous value, the monitor would automatically pump to a higher pressure by until it reaches a suitable pressure. Keep relax and wait for the measurement.	

SPECIFICATIONS

Model No. : TD-3250B/TD-3250C

Power source: Four 1.5V AA size alkaline batteries

Size of monitor w/o cuff: 137mm (L) x 90 mm (W) x 54 mm (H)

Memory: 352 measurement results with date& time

Power saving: Automatic power off if idle for 3 minutes

System operating condition: 10°C~40°C(50°F~104°F), below 85% R.H.

Monitor storage condition: -20°C~60°C(4°F ~ 140°F), below 95% R.H.

Strip storage condition: 4°C~40°C(39.2°F~104°F), below 85% R.H.

KETONE warning: glucose value over 240 mg/dL(13.3 mmol/L)

Measurement unit: mg/dL or mmol/L

Linear range: 20~600 mg/dL (1.1~33.3 mmol/L)

Precision: < 5 %(CV)

Accuracy: ± 15 mg/dL(0.83mmol/L) when glucose ≤ 75 mg/dL (4.17mmol/L);

$\pm 20\%$ when glucose ≥ 75 mg/dL(4.17mmol/L)

Pressure Range: 0-300 mmHg

Heart Rate Range: 40-199 beats per minute

Measurement unit: mmHg or KPa

Accuracy of Pressure: ± 3 mmHg or $\pm 2\%$ of reading

Accuracy of Heart rate: $\pm 4\%$ of reading

Maximum inflation pressure: 300 mmHg

The device has been certified to meet the electrical and safety requirements of: IEC 60601-1, EN 60601-1, IEC 61010-1, EN 61010-1, EN 61010-2-101, EN 60601-1-2, EN 61326

SYMBOLS INFORMATIO

Symbol	Referent
	<i>In vitro</i> diagnostic medical device
	Consult instructions for use
	Type BF equipment
	Temperature limitation
	Use by
	Batch code
	Manufacturer
SN	Serial number
	Caution, consult accompanying documents
	Authorized representative in the European Community
	Humidity limitation
CE ₀₁₂₃	CE mark
	collection for electrical and electronic equipment

FCC STATEMENT

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

- Consult the dealer or an experienced radio/TV technician for help.

Operation of this equipment is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

FCC RF Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

NOTE

The changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

At all other times, you should contact your health care professional for assistance.

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
Sumber: American Diabetes Association (2008). Standards of Medical Care in Diabetes. Diabetes Care, 31 (Supplement 1): S12–S54.

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Sumber: 1999 WHO/ISH guidelines for the management of hypertension
 *2: American Diabetes Association: The Diabetes-Heart Disease Link Surveying Attitudes, Knowledge and Risk (2002)

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
Sub-jenis TD-3250	Fungsi Penghantaran Data
<p data-bbox="75 227 223 258">TD-3250B</p>  <p>The diagram shows a handheld device (TD-3250B) on the left and a desktop computer (PC) on the right. A cable is connected between the device and the PC, indicating a wired data transfer connection.</p>	<p data-bbox="576 227 948 374">Anda boleh melihat keputusan pada KOMPUTER melalui sambungan kabel.</p>
<p data-bbox="75 579 223 609">TD-3250C</p>  <p>The diagram shows a handheld device (TD-3250C) on the left and a desktop computer (PC) on the right. The device has a signal icon above it, representing a wireless connection to the PC.</p>	<p data-bbox="576 579 959 763">Anda boleh melihat keputusan pada KOMPUTER melalui sambungan tanpa wayar yang menggunakan teknologi Bluetooth.</p>

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